

A white and gold handheld massage gun. The handle is white with a gold-colored band near the head. The head is white with a black base containing several massage heads. The handle features a digital display showing '12' and four control buttons: a red power button, a square button, a circular button, and a plus button.

Apply to DS-A15, DS-A16

Scraping is guided by the theory of meridians and Shu points (traditional Chinese medicine). It scrapes and rubs regularly on the body to make the skin smooth. There are red light candles or dark red bleeding points and the "four corners" of the changes, so as to promote blood circulation and penetrate Sha, regulate the qi and blood, dredging meridians and collaterals and improve the metabolism of the whole body and improve sub-health in an all-round way.

The instrument can suck up the fat layer of some parts of the body and stimulate the meridians and acupoints of the human body. At the same time, activate the cell vitality in the area, improve the thermal energy of the whole body, and enhance the deep fat through thermal energy. The temperature of fat makes it dissolve and metabolize, eliminates excess substances, and improves the tightness of the collagen fiber in the skin (epidermis). Dimensional curve, reduce fat and moisture, but also tighten the skin and shape the body.

- (1) those with bleeding or bleeding tendency; skin diseases and warts are prohibited
- (2) pregnancy, menstruation, skin damage, skin sensitivity, postoperative conditions, patients with high fever or critical diseases are prohibited
- (3) bleeding prone diseases, such as thrombocytopenia, leukemia, hemophilia, blood disorders
- (4) excessive fatigue, drunkenness and staying up too late are prohibited (2)

[illegible]

Product information

Model: _____

Color: _____

Purchase channel: _____

Customer information

Full name: _____

Address: _____


Date of purchase: _____

After inspection, the product meets the quality standards specified by the state and is delivered.

Product Name: Intelligent negative pressure cupping massage instrument

Product model: _____

Product color: _____

Quality testing:  QC
PASS
Check

Date of manufacture: _____

After starting the product, select the product function according to the needs, and start using it after applying essential oil to the required parts of the body. For sensitive skin, please test on the inner side of the arm or the root of the hair on the back of the hand before using it. If there is redness, swelling and itching, discontinue use immediately.

When scraping, apply some essential oil to the required parts (do not add more essential oil, let alone add essential oil to the machine). Turn the negative pressure pump switch on and start the negative pressure sucking the skin. The scraping and swabbing of each part is generally less than 20 minutes, with the strength from light to light/heavy, from small to large.

Note: the machine is equipped with large filter cotton, which can effectively prevent essential oil from entering the machine and plugging the protective filter cotton. When the filter cotton is dirty, please replace the filter cotton in time. The filter cotton can be recycled after cleaning and air drying. Negative pressure during cuppingTurn to the appropriate gear, suck the skin for 3-5 minutes, tilt the machine and pull it up. When the pressure relief key is easily pulled out and cuppingNo essential oil is required.

1. Please put the product in a safe place to avoid falling. Do not put heavy objects on the product to prevent damage.
2. Please store in a dry and well ventilated place to avoid strong sunlight, corrosive gas and liquid.
3. Please keep this product out of the reach of children.
4. When not in use for a long time, turn off the power supply and keep it fully charged every two months.
5. In case of technical failure, please contact technicians for maintenance in time. Do not disassemble the product without permission.
Before repair, the power supply must be cut off and the cleaning can be started only after the machine is cooled.
6. Clean with soft and slightly wet sponge or soft cloth. Do not clean the product with corrosive chemical cleaner.
7. This machine is not waterproof. It is strictly prohibited to use in humid environment, and the machine body must not be washed with water to avoid water or other liquids. The body enters the machine body to avoid short circuit or electric shock during use.

Technical parameter

Base

Control panel

Display window


Charging interface

Product Name: non-ferrous pressure-coupling stainless steel (type)
Material: Stainless steel ball, ABS, PC, etc. (steel ball)
Net weight: About 300g
Power: 8W
Charging parameters: 5V 2A
Rated voltage: 3.7V
Product size: 200 * 90 mm
Timing: 20 minutes
Note: Technical parameters are subject to change without notice

Timing: 20 minutes

Line after being fully charged, and cannot be started during charging
The white light flashes when charging and lights up when it's full. The digital
machine is equipped with a 20 minute timing function. The digital
display screen displays a 20 minute countdown and automatically
turns off the power after 20 minutes.
When the power switch is set up, press and hold for 2 seconds to
start the machine, the digital 00 on the digital display screen will
light up, and the 20 minute countdown will be displayed after 2
seconds; long press 2 seconds to stop it off, and the white light on the
digital display screen goes out. Short press the power key to stop

key game to restore the working state.
machine from working and enter the standby power.
(-) 1- gear strengthening / weaken keys, short press to adjust the
pressure, positive pressure. There are gears 0 - 21 gear, 0 is negative
pressure, 0 is the machine defaults 0 gears. Each time it is adjusted
(pressing) on the machine display shows the corresponding gears, and the
[pressure relief] it is effective when the negative pressure is turned
on. Press it firmly and unfold the gas for 2 seconds, which can easily
solve the problem of difficult turning.
(heating) short press to adjust the temperature. The heating function
(0.9 - 9.9 °C) is turned on, and the temperature is increased by 0.9 °C.
0.9 - 9.9 °C, 0.9 °C is the default value, the temperature is increased by
one degree, the red light of the classis lights up when the heating
function is turned on, and the temperature is about 55 degrees, which
is shorter when it reaches the 27th gear press to return to gear 0, and the
digital display screen will display the temperature gear, and the 20
minute countdown will be displayed again after 2 seconds, and the
classis has buttons friction relief, elastic steel balls can be used for hot
compress messages, 90 °C circulation and blood circulation, warming the
palm and protecting the stomach.

 This machine is stone needs charging type: it is fully charged without line operation, and cannot be started during charging (charging white light digital) 00 flashes and lights up after filling. The machine is equipped with a 20-minute automatic protection when working, and the digital display screen displays the 20 minute countdown, which will be automatically activated after 20 minutes stop working and turn off the power. If there is no operation within 5 minutes after startup, the machine will turn off the power automatically.

② When the power switch is set up, briefly press once to start the machine, the digital 00 on the digital display screen will light up, and the 5-minute countdown will be displayed after 5 seconds. If there is no operation, the machine will be automatically short press to shut down, and the white light on the digital display screen will go out.

[negative pressure] negative pressure gear adjustment key. Short press to adjust the negative pressure gear. It is the third gear in the default standby mode, there are gear 0-12 (gear 0 means no negative pressure).

[heating] short press to adjust the temperature. The heating function (i.e. gear 0) is not turned on by default when starting up. There are gears 0-12. Each gear is added a temperature increase by 0.5 degrees Celsius. When the heating function is turned on, the chassis red light is on, and the temperature is about 45 degrees [micro electricity] short press to adjust the micro current intensity.

The micro electricity function (i.e. gear 0) is not turned on by default when starting up. Press the micro electricity key on the steel ball at the bottom will generate micro current to stimulate human skin cells and tissues (not electric leakage). If there is any discomfort, stop of use).

[pressure relief] It is effective when the negative pressure is turned on. Press it briefly and hold the gas for 1 second, which can ease the problem of difficult tank lifting.

※ When the above negative pressure, heating and micro electric functions reach 12, press the corresponding function key again for a short time to turn off the function. When the digital display screen will display the corresponding functional gear, and the 20 minute countdown will be displayed again after 5 seconds.

※ When the above negative pressure, heating and micro electric functions reach gear 12, press the corresponding function key again for a short time to return to gear 0. The digital display will be displayed after each adjustment. The screen displays the corresponding functional gear, and the countdown of 20 minutes will be displayed again after 5 seconds.

咬口 蓝红黄黑 刮痧仪英文长款说明书天地翻 400x140mm 2022年7月22日 印前请仔细校对，印后恕不负责！